

LED Upgrade



 a|dec®

Upgrade to Brilliance

Bright, white, and true. Replacing your current A-dec light head with the new A-dec LED is quick, easy, and entirely hassle-free. Upgrade to 25% more illuminance while consuming one-fifth of the power needed for halogen bulbs.



Increase comfort, performance, and reduce your cost of ownership.

- See colors accurately for better tissue diagnosis
- Gain stunning white light with a color rendering index of 94 at 5,000K
- Reduce shadowing, improve contrast without strain to the eyes
- Optimize ergonomics for the entire dental team
- Position with three axes of rotation and 540° of horizontal motion
- Illuminate without curing with a cure-safe mode at 23,000 lux
- Lower energy costs by 80% and never replace a bulb again
- Treat patients in comfort—without radiating heat or noise

Am I eligible for an upgrade? As long as you have a working flexarm on an A-dec light manufactured after 1993, you can upgrade it to an A-dec LED.*

Will my new LED be compatible with my A-dec touchpad? Yes. A-dec touchpads utilizing DCS (Data Communication System) give you the choice to operate the light using the touchpad or local switching. With A-dec DCS, you can power on/off, select from three intensities, activate cure-safe mode, and support auto on/off with chair presets. Without DCS, your new A-dec LED light will have similar control via local switching only.

How can I learn more? Visit a-dec.com/LED or contact your Authorized A-dec Dealer.

Ordering Information

Description	Part Number
A-dec LED light head upgrade in White	28.1924.00
A-dec LED light head upgrade in Surf	28.1925.00

* Compatibility exceptions apply for the *A-dec Performer Light* and some *A-dec 300 lights*. For details, see your Authorized A-dec Dealer or contact an A-dec Customer Service Representative.



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132 USA
Tel: 1.800.547.1883 within USA/CAN
Tel: +1.503.538.7478 outside USA/CAN
Fax: +1.503.538.0276
a-dec.com

A-dec Australia

Tel: 1.800.225.010 within AUS
Tel: +61.(0)2.8332.4000 outside AUS

A-dec China

Tel: 400.600.5434 within CN
Tel: +86.571.89026088 outside CN

A-dec United Kingdom

Tel: 0800.233285 within UK
Tel: +44.(0).24.7635.0901 outside UK