

A-dec Lever Foot Control



“

Think of it as cruise control. Once you have determined the desired speed, you can relax your foot. You still have complete functionality, but drive in a more relaxed state.

Sam Zemke, Senior Product Manager

”

A-dec Lever Foot Control

A slight, natural shift of the foot delivers precise modulation of your electric handpiece motor speed. And once you reach the desired handpiece speed, you can rest your foot flat on the floor and continue working in a more relaxed state. A-dec's newest foot control is designed with straightforward left or right lever movement, which means you can intuitively switch between wet and dry cuts without taking your eyes off the patient.

Advantages:

- Precisely modulate and control electric handpiece motor speed.
- Rest foot during operation to relax foot and leg.
- Intuitive wet/dry side-to-side operation and chip blower/accessory button keeps focus on the patient.
- Encourages ergonomic “athletic stance” positioning.
- Compatible with A-dec 300 and 500 delivery systems.
- Backed by A-dec's 5-year warranty.

To learn more contact your local A-dec dealer, or visit a-dec.com.

